

MEDICAL TUBING CAPABILITIES



Balloons




Co-extrusions



Dual/Variable Flex



Microbore




Multi-layer




Multi-lumen




Over the Wire



Para-Tubing



Profile




Rod/Beading



Single-lumen



Striping



Tapered or Bumped

SECONDARY OPERATIONS

- Assembly
- Bonding
- Braiding or Braid Wire Reinforcement
- Coating
- Coiling or Coil Wire Reinforcement
- Cutting
- Flaring
- Insert Molding
- Laminating
- Laser Drilling
- Multi Durometer
- Multi Segment
- Overmolding
- Printing
- Punching
- Shaping
- Skiving
- Stretching
- Swaging
- Tipping

AVAILABLE MATERIALS

- Aesno (Nylon 12)
- Besno (Nylon 11)
- Braid - Aramid Fibers
- Braid - Round or SS Flat Wire
- Ethylene Vinyl Acetate (EVA)
- Fluoropolymers
- Fluoropolymer Composite (PTFE-Composite)
- Grilamid[®] (Nylon 12)
- Nanofiber Composite
- Pebax[®] (PEBA)
- Polyetheretherketone (PEEK)
- Pellethane[®] (Polyurethane)
- Polyester Terephthalate (PET)
- Polyamide (Nylon)
- Polycarbonate (PC)
- Polyethylene (PE)
- Polyimide (PI)
- Polyolefin (TPO)
- Polyphenylsulfone (PPSU)
- Polypropylene (PP)
- Polyurethane (PU)
- Polyacetal (POM)
- Polytetrafluoroethylene (PTFE)
- Polyurethane (PU, PUR, TP)
- Tecoflex[®] (Polyurethane)
- Texin[®] (TPE)
- Thermoplastic Elastomer
- Vestamid[®] (Nylon 12)

As the leading industry expert for a wide range of medical tubing, Vention offers a broad range of specialty tubing for the most unique and difficult applications. Vention Medical offers world class tolerances with industry leading process capabilities*:

- > 1.5C_{pk}
- i.d.'s as small as .004"
- walls as thin as .00015"
- o.d.'s as large as .5"

*Some capabilities are material dependent

Standard off the shelf tubing is available for delivery in as little as one day. Go to our online store at ventionmedical.com for over a thousand SKU's available for shipment today. Non-stock tubing can be delivered in as little as a few days to two weeks.

We can provide a full range of solutions from quick turn prototyping to high volume production. We will strategically partner with you to understand your needs and find the best solution. With over 30 years of experience, our quality is unparalleled. All facilities are ISO certified.

**CONTACT US TODAY FOR A QUOTE OR TO REQUEST
A SAMPLE AT INFO@VENTIONMEDICAL.COM**

ventionmedical.com